

# BIOCARRIER 650

2-8°C, 15-25°C



**SAMPLE  
TRANSPORT  
PACKAGING**



**The Intelsius BioCarrier range is designed for local courier routes and multiple site sample collection or drug protection.**

When used with 95kPa primaries, the BioCarrier 650 acts as a leak proof secondary and outer packaging in order to meet regulatory requirements. Intelsius recommends the use of 95kPa compliant pouches (available separately) if your primaries are not 95kPa.

## THIS PRODUCT CAN BE USED FOR

- Multisite sample collections
- Diagnostic samples
- Screening programmes

## FEATURES

- Qualified for both 2-8°C and 15-25°C this product offers thermal protection for both temperature ranges with the same bill of materials
- Supplied pre labelled for Category B, Biological Substances
- Rigid reusable transit case with leak proof EPP liner

## BIOCARRIER 650 RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions	Payload Volume	Mass	Volumetric weight L/6	2-8°C Duration (hours)	15-25°C Duration (hours)
BC022	BioCarrier 4 650 8.0L	397 x 295 x 326 mm 15.6 x 11.6 x 12.8 in	200 x 200 x 200 mm 7.8 x 7.8 x 7.8 in	8.0L	4.9 kg 10.8 lbs	6.36 kgvol <sup>c</sup> 14.0 lbsvol <sup>c</sup>	15:30	24:00
BC026	BioCarrier 5 650 10.8L	397 x 295 x 425 mm 15.6 x 11.6 x 16.7	200 x 200 x 270 mm 7.8 x 7.8 x 10.6 in	10.8L	6.6 kg 14.6 lbs	8.3 kgvol <sup>c</sup> 18.3 lbsvol <sup>c</sup>	19:50	24:00

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)

# BIOCARRIER 650

2-8°C & 15-25°C SOLUTIONS



**SAMPLE**  
TRANSPORT  
PACKAGING

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt Ltd  
2, Palak Pluto, Off Link Road,  
Kanderpada – Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 6525 6621  
f: +91 22 2892 3675  
compliance@intelsius.com