

# BIOSHIELD FK

VETERINARY FIELD KIT



HUMAN &  
ANIMAL  
HEALTH



The BioShield FK Veterinary Field Kit is a patented solution for the safe and effective transport of blood samples from extraction to laboratory presentation. The Category B compliant BioShield FK is equipped for sample traceability with sequential barcode numbering on all test tubes.

The BioShield FK is designed to give you peace of mind by providing your samples maximum protection during collection and transport. The inclusion of 95kPa primaries and leakproof secondaries ensures the kit is PI650 compliant.

## THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Screening programmes
- Eradication schemes

## FEATURES

- Sequentially numbered barcodes on all test tubes for sample traceability
- Complete Category B solutions including 95kPa test tubes, needles, absorbent and leak proof bag
- Kits available with a variety of needle gauges and lengths

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)

# BIOSHIELD FK

## VETERINARY FIELD KIT



HUMAN &  
ANIMAL  
HEALTH



### VETERINARY FIELD KIT RANGE

Code	Product	External Dimensions	Needle Size	Needles per Kit	Tube Size	Tubes per Kit	EDTA	Kits per Pack
FK011	Vet Field Kit 20g x 1" x 50 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	20g x 1"	50	6ml	50	No	10
FK012	Vet Field Kit 20g x 1.5" x 50 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	20g x 1.5"	50	6ml	50	No	10
FK013	Vet Field Kit 20g x 1" x 20 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	20g x 1"	20	6ml	20	No	10
FK014	Vet Field Kit 20g x 1.5" x 20 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	20g x 1.5"	20	6ml	20	No	10
FK015	Vet Field Kit (Tubes Only)	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	N/A	N/A	6ml	50	No	10
FK017	Vet Field Kit (Box Only)	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	N/A	N/A	N/A	N/A	N/A	10
FK018	Vet Field Kit 18g x 1" x 20 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	18g x 1"	20	6ml	20	No	10
FK019	Vet Field Kit 18g x 1" x 50 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	18g x 1"	50	6ml	50	No	10
FK020	Vet Field Kit 18g x 1.5" x 20 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	18g x 1.5"	20	6ml	20	No	10
FK021	Vet Field Kit 18g x 1.5" x 50 x 10	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	18g x 1.5"	50	6ml	50	No	10
FK022	Vet Field Kit (Tubes Only)	310 x 150 x 135 mm 12.2 x 5.9 x 5.3 in	N/A	N/A	6ml	20	No	10

Product Name Key:[Range] [Needle Weight x Needle Length x Tube Quantity x Kit Quantity]

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt Ltd  
2, Palak Pluto, Off Link Road,  
Kanderpada - Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 6525 6621  
f: +91 22 2892 3675  
compliance@intelsius.com