



The Intelsius BioShield Range of unique patented systems are designed to meet all handling processes from sample extraction to laboratory presentation. BioShields are highly protective single-use shipping solutions. They are available for a variety of sample capacities, helping to reduce waste and cut costs associated with sample transportation.

The systems are well suited to a wide variety of Category B samples, including solids, liquids, and samples held in fixative. The BioShield helps reduce total costs by ensuring that all samples will arrive at the same time for a simple and efficient sample reception.

## THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Screening programmes
- Eradication schemes

## FEATURES

- Complete Category B compliant solutions including 95kPa 90ml sample pots, tray, extraction utensils and leak proof bag
- Techni-Ice gel sheet available with next day systems for temperature stability during transport

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



## BIOSHIELD 'SAME DAY' RANGE

| Product         | Product Code (Inc. Pots) | Product Code (Exc. Pots) | Product Code (Inc. Pots, Spoon & Forceps) | Sample Capacity | External Dimensions                        | Volumetric Weight                                  |
|-----------------|--------------------------|--------------------------|---|-----------------|--|--|
| BioShield 6     | BSE061                   | BSE072                   | BSE055                                    | 6               | 291 x 191 x 102 mm<br>11.5 x 7.6 x 4.1 in  | 0.94 kgvol <sup>c</sup><br>2.1 lbsvol <sup>c</sup> |
| BioShield 10    | BSE075                   | BSE077                   | BSE038                                    | 10              | 291 x 191 x 102 mm<br>11.5 x 7.6 x 4.1 in  | 0.94 kgvol <sup>c</sup><br>2.1 lbsvol <sup>c</sup> |
| BioShield SR 20 | BSE056                   | BSE065                   | BSE057                                    | 20              | 519 x 340 x 125 mm<br>20.5 x 13.4 x 5 in   | 3.7 kgvol <sup>c</sup><br>8.1 lbsvol <sup>c</sup>  |
| BioShield SR 40 | BSE058                   | BSE066                   | BSE059                                    | 40              | 519 x 340 x 125 mm<br>20.5 x 13.4 x 5 in   | 3.7 kgvol <sup>c</sup><br>8.1 lbsvol <sup>c</sup>  |
| BioShield TR 80 | BSE083                   | BSE021                   | BSE081                                    | 80              | 495 x 355 x 245 mm<br>19.5 x 14.0 x 9.7 in | 7.2 kgvol <sup>c</sup><br>15.9 lbsvol <sup>c</sup> |

## BIOSHIELD 'NEXT DAY' RANGE

| Product         | Product Code (Inc. Pots) | Product Code (Exc. Pots) | Product Code (Inc. Pots, Spoon & Forceps) | Sample Capacity | External Dimensions                        | Volumetric Weight                                  |
|-----------------|--------------------------|--------------------------|---|-----------------|--|--|
| BioShield 6     | BSE062                   | BSE071                   | BSE051                                    | 6               | 291 x 191 x 102 mm<br>11.5 x 7.6 x 4.1 in  | 0.94 kgvol <sup>c</sup><br>2.1 lbsvol <sup>c</sup> |
| BioShield 10    | BSE074                   | BSE076                   | BSE073                                    | 10              | 291 x 191 x 102 mm<br>11.5 x 7.6 x 4.1 in  | 0.94 kgvol <sup>c</sup><br>2.1 lbsvol <sup>c</sup> |
| BioShield SR 20 | BSE050                   | BSE063                   | BSE064                                    | 20              | 519 x 340 x 125 mm<br>20.5 x 13.4 x 5 in   | 3.7 kgvol <sup>c</sup><br>8.1 lbsvol <sup>c</sup>  |
| BioShield SR 40 | BSE041                   | BSE042                   | BSE067                                    | 40              | 519 x 340 x 125 mm<br>20.5 x 13.4 x 5 in   | 3.7 kgvol <sup>c</sup><br>8.1 lbsvol <sup>c</sup>  |
| BioShield TR 80 | BSE019                   | BSE007                   | BSE082                                    | 80              | 495 x 355 x 245 mm<br>19.5 x 14.0 x 9.7 in | 7.2 kgvol <sup>c</sup><br>15.9 lbsvol <sup>c</sup> |

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt Ltd  
2, Palak Pluto, Off Link Road,  
Kanderpada - Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 6525 6621  
f: +91 22 2892 3675  
compliance@intelsius.com