

BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



The BioTherm range of Category A and Category B solutions come complete with the Intelsius PathoPak sample transport systems and are fully compliant to IATA, ADR and DOT (49 CFR) transport regulations. The systems are qualified to ensure sample integrity for up to 96 hours against 2-8°C, 15-25°C and dry ice. Each system is provided with handles and the correct compliance labelling to enable careful handling throughout the cold chain.

THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Complete temperature controlled Category A and Category B compliant overpack systems
- Range of systems to transport samples at 2-8°C, 15-25°C and deep frozen to maintain sample integrity
- 4GU compliant packaging systems for use with any primary receptacle
- Ergonomically engineered, reusable, sample bottle for safe manual handling and easy extraction

For more information please contact Intelsius UK cs@intelsius.com

BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



DRY ICE SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Empty Mass	Dry Ice Mass	Volumetric Weight L/6	Duration
BT090	BioTherm 10A Dry Ice	357 x 309 x 315 mm 14.1 x 12.2 x 12.5 in	90 dia x 170 mm 3.5 dia x 6.7 in	1 L	1.4 kg 3.1 lbs	6.9 kg 15.3 lbs	5.8 kgvol ^c 12.8 lbsvol ^c	72
BT052	BioTherm 15A Dry Ice	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	2.4 kg 5.3 lbs	11.5 kg 25.4 lbs	11.1 kgvol ^c 24.4 lbsvol ^c	96
BT062	BioTherm 30A Dry Ice	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 255 mm 4.9 dia x 10.0 in	3 L	3.3 kg 7.1 lbs	25.0 kg 55.2 lbs	17.0 kgvol ^c 37.4 lbsvol ^c	96 +

HIGH PERFORMANCE 2-8°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT056P	BioTherm 30A HP 2-8 96 Summer	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	10.2 kg 22.4 lbs	17 kgvol ^c 37.4 lbsvol ^c	96
BT057P	BioTherm 30A HP 2-8 96 Winter	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	9.0 kg 19.8 lbs	17 kgvol ^c 37.4 lbsvol ^c	96

CRT 15-25°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT105	BioTherm 30A CRT 15-25	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	11.9 kg 26.1 lbs	17 kgvol ^c 37.4 lbsvol ^c	96

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt Ltd
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 6525 6621
f: +91 22 2892 3675
compliance@intelsius.com