

BIO THERM EXTREME 6

1-10°C



**SAMPLE
TRANSPORT
PACKAGING**



Optimised to deliver 36 hours of thermal protection, the BioTherm Extreme 6 is the latest addition to the Intelsius BioTherm Extreme range. This easy to use water, EPS & VIP hybrid system is ideal for UK, Europe and North America domestic and overnight shipments and ensures that your biological samples are protected against extreme climates.

FEATURES

- 36 hours performance qualified against ISTA 7D HOT 1-10°C
- Hybrid system of water, EPS & VIP
- Optimised for performance vs shipping costs: Sub 3kg volumetric weight and sub 3kg physical weight
- Payload space for one cryovial box or up to 12 blood tubes
- Category B compliant system when used in conjunction with PathoSeal 95 (available separately)
- Ideal for UK, Europe and North America domestic/overnight shipments
- Easy to assemble, utilising one way fit coolant cassette

Product Code	Product Name	External Dimensions	Payload Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
EES6	BioTherm Extreme 6	280 x 275 x 180 mm 11.0 x 10.8 x 7.1 in	165 x 165 x 62 mm 6.5 x 6.5 x 2.4 in	1.6 L	2.4 kg 5.3 lbs	2.8 kgvol ^c 6.1 lbsvol ^c	36

For more information please contact Intelsius UK cs@intelsius.com

BIO THERM EXTREME 6

1-10°C



SAMPLE
TRANSPORT
PACKAGING

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt Ltd
2, Palak Pluto, Off Link Road,
Kanderpada – Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 6525 6621
f: +91 22 2892 3675
compliance@intelsius.com