

PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



The PathoPak range provides strong, durable and reusable, rigid secondary containers, compliant with transport regulations for UN Class 6.2 Infectious Substances Category A and B. Available in three sizes 1, 2 and 3 Litres.

THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES:

- Complete Category A and Category B solutions include 95kPa tested secondary container, absorbent, bubble wrap, labels, and UN markings
- Reusable
- Ergonomically engineered bottle design for safe manual handling including 1.5 turn lid and non-tapered aperture for easy extraction of samples
- 4GU certified for the use of any primary
- Suitable for cleaning by autoclave
- Full refurbishment services available

For more information please contact Intelsius UK cs@intelsius.com



PATHOPAK RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP011	DGP PathoPak 1L	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol ^c 0.86 lbsvol ^c
PP012	DGP PathoPak 2L	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol ^c 1.69 lbsvol ^c
PP013	DGP PathoPak 3L	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol ^c 2.63 lbsvol ^c

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt Ltd
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 6525 6621
f: +91 22 2892 3675
compliance@intelsius.com