

PATHOSHIELD



**SAMPLE
TRANSPORT
PACKAGING**



The PathoShield range are complete shipping solutions designed to carry a wide range of primaries, including blood vials, swabs and sample containers. Each system is pre-printed with the required markings and includes the relevant components to ship samples in compliance with transport regulations.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Complete Category A and Category B solutions include flexible 95kPa tested compliant secondary container, absorbent, labels, and UN markings
- Range of sizes suitable for shipping through postal and courier networks
- Carry a wide range of primaries including blood vials, swabs and sample containers
- Can be supplied as an insulated pack
- Supplied flat packed for efficient storage and shipping prior to use

For more information please contact Intelsius UK cs@intelsius.com



PATHOSHIELD RANGE

Product Code	Product Name	External Dimensions	Pouch Internal Dimensions (W x L)	Number of Pouches	Number of Tubes	Mass	Volumetric Weight L/6
--------------	--------------	---------------------	-----------------------------------	-------------------	-----------------	------	-----------------------

Category A Solutions

PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	3	0.23 kg 0.5 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
-------	--------------------------	---	------------------------------	---	---	--------------------	---

Category B Solutions

PA003	PathoShield 3 w A4 Pouch	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	6	0.16 kg 0.35 lbs	0.5 kgvol ^c 1.0 lbsvol ^c
PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	228 x 250 mm 8.9 x 9.8 in	1	6	0.23 kg 0.5 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA001	PathoShield 5 w A5 Pouch	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.09 kg 0.20 lbs	0.4 kgvol ^c 0.8 lbsvol ^c
PA006	PathoShield 6 w A5 Pouch	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.07 kg 0.15 lbs	0.2 kgvol ^c 0.4 lbsvol ^c
PA007	PathoShield 7 w A5 Pouch	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	155 x 180mm 6.1 x 7.0 in	1	4	0.10 kg 0.21 lbs	0.3 kgvol ^c 0.6 lbsvol ^c
PA008	PathoShield 8 w A5 Pouch	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	155 x 180mm 6.1 x 7.0 in	1	4	0.04 kg 0.09 lbs	0.1 kgvol ^c 0.2 lbsvol ^c
PA023	PathoShield 3 w A4 Bag	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.13 kg 0.28 lbs	0.5 kgvol ^c 1.0 lbsvol ^c
PA024	PathoShield 4 w A4 Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.19 kg 0.43 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA025	PathoShield 5 w A5 Bag	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.08 kg 0.18 lbs	0.4 kgvol ^c 0.8 lbsvol ^c
PA026	PathoShield 6 w A5 Bag	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.06 kg 0.13 lbs	0.2 kgvol ^c 0.4 lbsvol ^c
PA027	PathoShield 7 w A5 Bag	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.09 kg 0.19 lbs	0.3 kgvol ^c 0.6 lbsvol ^c
PA028	PathoShield 8 w A5 Bag	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.04 kg 0.09 lbs	0.1 kgvol ^c 0.2 lbsvol ^c



PATHOSHIELD 4 INSULATED RANGE

Product Code	Product Name	External Dimensions	Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6
BT021	PathoShield 4 Insulated	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	251 x 151 x 62 mm 9.8 x 5.9 x 2.4 in	2.3 L	0.24 kg 0.53 lbs	0.9 kgvol ^c 2.1 lbsvol ^c

Product Code	Product Name	External Dimensions	Internal Dimensions	Number of Pouches	Number of Tubes	Mass	Volumetric Weight L/6
PA037	PathoShield 4 Insulated w A5 Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	150 x 190 mm 5.9 x 7.4 in	1	6	0.26 kg 0.57 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA038	PathoShield 4 Insulated w A5 Bag & Gel Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	150 x 190 mm 5.9 x 7.4 in	1	6	0.66 kg 1.45 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA039	PathoShield 4 Insulated w A5 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	155 x 180 mm 6.1 x 7.0 in	1	6	0.26 kg 0.57 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA040	PathoShield 4 Insulated w A5 Pouch & Gel Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	155 x 180 mm 6.1 x 7.0 in	1	6	0.66 kg 1.45 lbs	0.9 kgvol ^c 2.1 lbsvol ^c

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt Ltd
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 6525 6621
f: +91 22 2892 3675
compliance@intelsius.com