

PATHOSORB ABSORBENT



COMPONENT
PACKAGING



The Intelsius PathoSorb is a Super Absorbent Polymer (SAP) that ensures excellent absorption capacity to size ratio. Upon absorption during accidental spillage, PathoSorb acquires a 'gel like' consistency, locking in the spilled material for easy disposal. PathoSorb is lint free, reducing sterile contamination risk.

Category A and Category B transport regulations specify that absorbent material must be placed between primary and secondary containers in sufficient quantity to absorb the entire contents of the primary containers.

FEATURES

- PathoSorb utilises SAP technology to provide effective absorbency whilst maximising available payload space.
- For multiple sample transportation, PathoSorb 4 & 6 Bay pouches provide a dual purpose; simultaneously protecting and separating samples.

For more information please contact Intelsius UK cs@intelsius.com



PATHOSORB ABSORBENT



COMPONENT
PACKAGING



PATHOSORB ABSORBENT RANGE

Product Code	Product Name	Dimensions	Water Absorbent Capacity	Saline Absorbent Capacity
AB001	PathoSorb 50ml	150 x 75 mm 6 x 3 in	56 ml 1.9 fl oz	50 ml 1.7 fl oz
AB010	PathoSorb 100ml	160 x 140 mm 6.3 x 5.5 in	112 ml 3.9 fl oz	100 ml 3.5 fl oz
AB009	PathoSorb 200ml	225 x 175 mm 8.9 x 6.9 in	223 ml 7.8 fl oz	200 ml 7 fl oz
AB002	PathoSorb 300ml	305 x 220 mm 12 x 8.7 in	335 ml 11.7 fl oz	300 ml 10.5 fl oz
AB007	PathoSorb 600ml	470 x 295 mm 18.5 x 11.6 in	693 ml 24.3 fl oz	600 ml 21 fl oz
AB005	PathoSorb 4 Bay Pouch	225 x 125 mm 8.9 x 4.9 in	169 ml 5.9 fl oz	152 ml 5.3 fl oz
AB008	PathoSorb 6 Bay Pouch	325 x 125 mm 12.8 x 4.9 in	283 ml 8.3 fl oz	215 ml 7.5 fl oz

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt Ltd
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 6525 6621
f: +91 22 2892 3675
compliance@intelsius.com