

# PATHOSEAL

## SAMPLE TRANSPORT BAG



COMPONENT  
PACKAGING



**Intelsius offer a range of flexible secondary packaging ideal for the transportation of a variety of primary containers such as blood vials, swabs, blood collection bags and medical devices.**

The PathoSeal 95 range is designed and tested to withstand a pressure difference of 95 kPa to meet UN3373 Biological Substance, Category B regulations.

The PathoSeal Liquitite range is designed and tested to meet Category B transport regulations when used in conjunction with 95 kPa primaries.

All PathoSeal products feature a simple closure system enabling ease of packing, plus the unique tamper evident seal ensures the integrity of samples. In addition the bags are sequentially numbered, providing a chain of custody and additional security. An external document pouch is also present for all sizes with the exception of PathoSeal 95 A6.

### THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Category B, Biological Substances
- Clinical trials
- Screening programmes
- Eradication schemes

### FEATURES

- Suitable for use within Category B solutions
- Tamper evident seal
- Sequentially numbered for sample traceability

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



# PATHOSEAL

## SAMPLE TRANSPORT BAG



COMPONENT  
PACKAGING



### PATHOSEAL 95 RANGE

Product Code	Product Name	Internal Dimensions (L x W)
BZ014	PathoSeal 95 A3	394 x 310 mm 15.5 x 12.2 in
BZ015	PathoSeal 95 A3 with 300ml Absorbent	394 x 310 mm 15.5 x 12.2 in
BZ006	PathoSeal 95 A4	284 x 200 mm 11.2 x 7.9 in
BZ009	PathoSeal 95 A4 with 200ml Absorbent	284 x 200 mm 11.2 x 7.9 in
BZ007	PathoSeal 95 A5	201 x 152 mm 7.9 x 6.0 in
BZ008	PathoSeal 95 A5 with 50ml Absorbent	201 x 152 mm 7.9 x 6.0 in
BZ020	PathoSeal 95 A6	139 x 99 mm 5.5 x 3.9 in
BZ013	PathoSeal 95 A6 with 50ml Absorbent	139 x 99 mm 5.5 x 3.9 in

### PATHOSEAL LIQUITITE

Product Code	Product Name	Internal Dimensions (L x W)
BZ001	PathoSeal Liquitite 29/18	720 x 450 mm 28.3 x 17.7 in

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt Ltd  
2, Palak Pluto, Off Link Road,  
Kanderpada - Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 6525 6621  
f: +91 22 2892 3675  
compliance@intelsius.com