

BIO THERM EXTREME 6

1-10°C



SAMPLE
TRANSPORT
PACKAGING



Optimised to deliver 36 hours of thermal protection, the BioTherm Extreme 6 is the latest addition to the Intelsius BioTherm Extreme range. This easy to use water, EPS & VIP hybrid system is ideal for UK, Europe and North America domestic and overnight shipments and ensures that your biological samples are protected against extreme climates.

FEATURES

- 36 hours performance qualified against ISTA 7D HOT 1-10°C
- Hybrid system of water, EPS & VIP
- Optimised for performance vs shipping costs: Sub 3kg volumetric weight and sub 3kg physical weight
- Payload space for one cryovial box or up to 12 blood tubes
- Category B compliant system when used in conjunction with PathoSeal 95 (available separately)
- Ideal for UK, Europe and North America domestic/overnight shipments
- Easy to assemble, utilising one way fit coolant cassette

| Product Code | Product Name | External Dimensions | Payload Dimensions | Payload Volume | Mass | Volumetric Weight L/6 | Duration |
|--------------|--------------------|--|---|----------------|-------------------|---|----------|
| EES 6 | BioTherm Extreme 6 | 280 x 275 x 180 mm 11.0 x 10.8 x 7.1 in | 165 x 165 x 62 mm 6.5 x 6.5 x 2.4 in | 1.6 L | 2.4 kg 5.3 lbs | 2.8 kgvol ^c 6.1 lbsvol ^c | 36 |

For more information please contact Intelsius UK cs@intelsius.com



BIO THERM EXTREME 6

1-10°C



SAMPLE
TRANSPORT
PACKAGING

DGP Intelsius LLC

7696 Zionsville Road,
Indianapolis
IN 46268
United States of America

AMERICAS

t: +1 317 452 4006
f: +1 317 290 4621

americasorders@intelsius.com

DGP Intelsius Ltd

1 Harrier Court
Airfield Business Park
Elvington
York, YO41 4EA
United Kingdom

EUROPE

t: +44 1904 607 390
f: +44 1904 607 391

cs@intelsius.com

DGP Intelsius Cold Chain Pvt Ltd

2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West)
Mumbai 400 068
Maharashtra, India

INDIA

t: +91 22 6525 6621
f: +91 22 2892 3675

compliance@intelsius.com

DGP Intelsius Sdn Bhd

Plot 11, Lorong Bemban 2
Estet Perindustrian Bemban
31000 Batu Gajah
Malaysia

ASIA

t: +605 3666 399
f: +605 3669 399

compliance@intelsius.com