

PATHOPOUCH 95

SAMPLE TRANSPORT POUCH



COMPONENT
PACKAGING



The PathoPouch 95 flexible secondary range is designed and tested to withstand a pressure difference of 95 kPa to meet Category A and Category B transport regulations. They are leak-proof and suitable to hold a wide range of primary containers, such as blood vials, sample containers, swabs, blood collection bags and medical devices.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Category A and B Infectious Substances
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Suitable for use within Category A and Category B solutions
- Tamper evident seal
- External document pouch
- Carry a wide range of primaries
- Available pre-packed with absorbent

For more information please contact Intelsius UK cs@intelsius.com

PATHOPOUCH 95

SAMPLE TRANSPORT POUCH



COMPONENT
PACKAGING



PATHOPOUCH 95 RANGE

Product Code	Product Name	Internal Dimensions (L x W)	Absorbent
BP007	PathoPouch Xcel	460 x 334 mm 18.1 x 13.1 in	-
BP008	PathoPouch Xcel with Absorbent	460 x 334 mm 18.1 x 13.1 in	6 Bay
BP005	PathoPouch A3	359 x 261 mm 14.1 x 10.3 in	-
BP006	PathoPouch A3 with Absorbent	359 x 261 mm 14.1 x 10.3 in	6 Bay
BP003	PathoPouch A4	308 x 230 mm 12.1 x 9.1 in	-
BP001	PathoPouch A4 with Absorbent	308 x 230 mm 12.1 x 9.1 in	6 Bay
BP004	PathoPouch A5	230 x 152 mm 9.1 x 6.0 in	-
BP002	PathoPouch A5 with Absorbent	230 x 152 mm 9.1 x 6.0 in	4 Bay

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom

t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America

t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia

t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland

t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India

t: +91 22 2892 3675
f: +91 22 2892 3675
compliance@intelsius.com