

PATHOSORB ABSORBENT



COMPONENT
PACKAGING



The Intelsius PathoSorb is a Super Absorbent Polymer (SAP) that ensures excellent absorption capacity to size ratio. Upon absorption during accidental spillage, PathoSorb acquires a 'gel like' consistency, locking in the spilled material for easy disposal. PathoSorb is lint free, reducing sterile contamination risk.

Category A and Category B transport regulations specify that absorbent material must be placed between primary and secondary containers in sufficient quantity to absorb the entire contents of the primary containers.

FEATURES

- ▶ PathoSorb utilises SAP technology to provide effective absorbency whilst maximising available payload space.
- ▶ For multiple sample transportation, PathoSorb 4 & 6 Bay pouches provide a dual purpose; simultaneously protecting and separating samples.

For more information please contact Intelsius UK cs@intelsius.com



PATHOSORB ABSORBENT



COMPONENT
PACKAGING



PATHOSORB ABSORBENT RANGE

| Product Code | Product Name | Dimensions | Water Absorbent Capacity | Saline Absorbent Capacity |
|--------------|--------------------------|--------------------------------|--------------------------|---------------------------|
| AB001 | PathoSorb 50ml | 150 x 75 mm 6 x 3 in | 56 ml 1.9 fl oz | 50 ml 1.7 fl oz |
| AB010 | PathoSorb 100ml | 160 x 140 mm 6.3 x 5.5 in | 112 ml 3.9 fl oz | 100 ml 3.5 fl oz |
| AB009 | PathoSorb 200ml | 225 x 175 mm 8.9 x 6.9 in | 223 ml 7.8 fl oz | 200 ml 7 fl oz |
| AB002 | PathoSorb 300ml | 305 x 220 mm 12 x 8.7 in | 335 ml 11.7 fl oz | 300 ml 10.5 fl oz |
| AB007 | PathoSorb 600ml | 470 x 295 mm 18.5 x 11.6 in | 693 ml 24.3 fl oz | 600 ml 21 fl oz |
| AB005 | PathoSorb 4 Bay Pouch | 225 x 125 mm 8.9 x 4.9 in | 169 ml 5.9 fl oz | 152 ml 5.3 fl oz |
| AB008 | PathoSorb 6 Bay Pouch | 325 x 125 mm 12.8 x 4.9 in | 283 ml 8.3 fl oz | 215 ml 7.5 fl oz |

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 2892 3675
f: +91 22 2892 3675
compliance@intelsius.com