

SAMPLE EXTRACTION UTENSILS





COMPONENT
PACKAGING



Ergonomically designed, each Sample Extraction Utensil is available in either straight or curved form, depending on your preference. The disposable tools are made with flexible and durable polypropylene material. These sample extraction utensils are exclusively designed by Intelsius.

THIS PRODUCT CAN BE USED FOR

-  BSE and Scrapie sample extraction
-  Extraction of soft tissue with minimal damage

For more information please contact Intelsius UK cs@intelsius.com








SAMPLE EXTRACTION UTENSILS



COMPONENT
PACKAGING

SAMPLE EXTRACTION UTENSILS (SEU) RANGE

Product Image	Code	Product
	SP001	BSE Curved SEU
	SP003	Scrapie Spatula SEU
	SP008	BSE Straight SEU
	SP009	Ovine Scrapie Straight SEU
	SP012	Forceps SEU

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 2892 3675
f: +91 22 2892 3675
compliance@intelsius.com