

# BIOMAILER



**SAMPLE  
TRANSPORT  
PACKAGING**



The Intelsius BioMailer is a Category B compliant transport solution designed for use within national and global postal networks to ship small samples such as blood vials and slides.

## THIS PRODUCT CAN BE USED FOR

- Category B, Biological Substances
- Screening programmes
- Clinical trials
- Swabs
- Vials

## FEATURES

- Reusable rigid plastic case ensures samples are protected during transit
- Refurbishment kits available to replace single use components

## BIOMAILER RANGE

Product Code	Product Name	External Dimensions	Internal Pouch Dimensions	Mass	Volumetric Weight L/6	Case Quantity
BM003	BioMailer with A5 Pouch	215 x 165 x 30 mm 8.5 x 6.5 x 1.2 in	155 x 180mm 6.1 x 7.0 in	0.2 kg 0.44 lbs	0.2 kgvol <sup>c</sup> 0.44 lbsvol <sup>c</sup>	50
BM004	BioMailer A5 Pouch Refurbishment Kit	N/A	155 x 180mm 6.1 x 7.0 in	0.04 kg 1.1 lbs	N/A	100
BM001	BioMailer with A5 Bag	215 x 165 x 30 mm 8.5 x 6.5 x 1.2 in	240 x 145mm 5.7 x 9.4 in	1.2 kg 0.44 lbs	0.2 kgvol <sup>c</sup> 0.44 lbsvol <sup>c</sup>	50
BM002	BioMailer A5 Bag Refurbishment Kit	N/A	240 x 145mm 5.7 x 9.4 in	0.04 kg 0.1 lbs	N/A	100

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)

# BIOMAILER



**SAMPLE**  
TRANSPORT  
PACKAGING

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
2, Palak Pluto, Off Link Road,  
Kanderpada - Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 2892 3675  
f: +91 22 2892 3675  
compliance@intelsius.com