



The Intelsius BioShield Range of unique patented systems are designed to meet all handling processes from sample extraction to laboratory presentation. BioShields are highly protective single-use shipping solutions. They are available for a variety of sample capacities, helping to reduce waste and cut costs associated with sample transportation.

The systems are well suited to a wide variety of Category B samples, including solids, liquids, and samples held in fixative. The BioShield helps reduce total costs by ensuring that all samples will arrive at the same time for a simple and efficient sample reception.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Screening programmes
- Eradication schemes

FEATURES

- Complete Category B compliant solutions including 95kPa 90ml sample pots, tray, extraction utensils and leak proof bag
- Techni-Ice gel sheet available with next day systems for temperature stability during transport

For more information please contact Intelsius UK cs@intelsius.com



BIOSHIELD 'SAME DAY' RANGE

Product	Product Code (Inc. Pots)	Product Code (Exc. Pots)	Product Code (Inc. Pots, Spoon & Forceps)	Sample Capacity	External Dimensions	Volumetric Weight
BioShield 6	BSE061	BSE072	BSE055	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield 10	BSE075	BSE077	BSE038	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield SR 20	BSE056	BSE065	BSE057	20	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield SR 40	BSE058	BSE066	BSE059	40	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield TR 80	BSE083	BSE021	BSE081	80	495 x 355 x 245 mm 19.5 x 14.0 x 9.7 in	7.2 kgvol ^c 15.9 lbsvol ^c

BIOSHIELD 'NEXT DAY' RANGE

Product	Product Code (Inc. Pots)	Product Code (Exc. Pots)	Product Code (Inc. Pots, Spoon & Forceps)	Sample Capacity	External Dimensions	Volumetric Weight
BioShield 6	BSE062	BSE071	BSE051	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield 10	BSE074	BSE076	BSE073	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield SR 20	BSE050	BSE063	BSE064	20	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield SR 40	BSE041	BSE042	BSE067	40	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield TR 80	BSE019	BSE007	BSE082	80	495 x 355 x 245 mm 19.5 x 14.0 x 9.7 in	7.2 kgvol ^c 15.9 lbsvol ^c

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom

t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America

t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia

t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland

t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
2, Palak Pluto, Off Link Road,
Kanderpada - Dahisar (West),
Mumbai 400 068, Maharashtra, India

t: +91 22 2892 3675
f: +91 22 2892 3675
compliance@intelsius.com