

# PATHOPAK



**SAMPLE  
TRANSPORT  
PACKAGING**



The PathoPak range provides strong, durable and reusable, rigid secondary containers, compliant with transport regulations for UN Class 6.2 Infectious Substances Category A and B. Available in three sizes 1, 2 and 3 Litres.

## THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

## FEATURES:

- Complete Category A and Category B solutions include 95kPa tested secondary container, absorbent, bubble wrap, labels, and UN markings
- Reusable
- Ergonomically engineered bottle design for safe manual handling including 1.5 turn lid and non-tapered aperture for easy extraction of samples
- 4GU certified for the use of any primary
- Suitable for cleaning by autoclave
- Full refurbishment services available

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)

# PATHOPAK



**SAMPLE  
TRANSPORT  
PACKAGING**



## PATHOPAK RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP011	DGP PathoPak 1L	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol <sup>c</sup> 0.86 lbsvol <sup>c</sup>
PP012	DGP PathoPak 2L	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol <sup>c</sup> 1.69 lbsvol <sup>c</sup>
PP013	DGP PathoPak 3L	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol <sup>c</sup> 2.63 lbsvol <sup>c</sup>

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com



**IRELAND:** DGP Intelsius Teo  
Unit 8, Ballyvourney Industrial Estate,  
Ballyvourney, Co. Cork,  
P12 A527, Ireland  
t: +353 (0) 26 456 99  
f: +353 (0) 26 456 98  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
2, Palak Pluto, Off Link Road,  
Kanderpada - Dahisar (West),  
Mumbai 400 068, Maharashtra, India  
t: +91 22 2892 3675  
f: +91 22 2892 3675  
compliance@intelsius.com