

# COVID-19 SAMPLE TRANSPORT PACKAGING SOLUTIONS



**SAMPLE  
TRANSPORT  
PACKAGING**

## COVID-19 PACKAGING

As a global supplier of compliant and certified sample transport packaging solutions, we have ensured placing orders and using our UN3373 Category B packaging is as simple as possible. All DGP Intelsius UN3373 solutions are supplied ready to ship COVID-19 samples for testing.

Intelsius are well positioned to cope with any sudden increase in demand due to the possibility of a global pandemic. Stock levels of all our UN3373 Category B packaging have been increased at all global sites to ensure demand can be satisfied and potentially infectious samples can be transported safely and compliantly.



## PATHOPAK

The PathoPak range provides strong, durable and reusable, rigid secondary containers, compliant with transport regulations for UN Class 6.2 Infectious Substances Category A and B. Available in three sizes 1, 2 and 3 Litres. The PathoPak is the recommended solution for shipping nasopharyngeal and oropharyngeal swabs.

- Complete Category A and Category B solutions include 95kPa tested secondary container, absorbent, bubble wrap, labels, and UN markings
- Ergonomically engineered bottle design for safe manual handling for easy extraction of samples



Product Code	Product Name	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP011	DGP PathoPak 1L	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol <sup>c</sup> 0.86 lbsvol <sup>c</sup>
PP012	DGP PathoPak 2L	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol <sup>c</sup> 1.69 lbsvol <sup>c</sup>
PP013	DGP PathoPak 3L	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol <sup>c</sup> 2.63 lbsvol <sup>c</sup>



For more information please contact Intelsius UK [compliance@intelsius.com](mailto:compliance@intelsius.com)

# COVID-19 SAMPLE TRANSPORT PACKAGING SOLUTIONS



**SAMPLE  
TRANSPORT  
PACKAGING**

## PATHOSHIELD

The PathoShield range are complete shipping solutions designed to carry a wide range of primaries, including blood vials, swabs and sample containers. Each system is pre-printed with the required markings and includes the relevant components to ship samples in compliance with transport regulations.

- Range of sizes suitable for shipping through postal and courier networks
- Carry a wide range of primaries including blood vials, swabs and sample containers

Product Code	Product Name	External Dimensions	Pouch Internal Dimensions (W x L)	Mass	Volumetric Weight L/6
<b>Category B Solutions</b>					
PA028	PathoShield 8 w A5 Bag	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	150 x 190 mm 5.9 x 7.4 in	0.04 kg 0.09 lbs	0.1 kgvol <sup>l</sup> 0.2 lbsvol <sup>l</sup>
PA026	PathoShield 6 w A5 Bag	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	0.06 kg 0.13 lbs	0.2 kgvol <sup>l</sup> 0.4 lbsvol <sup>l</sup>
PA027	PathoShield 7 w A5 Bag	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	150 x 190 mm 5.9 x 7.4 in	0.09 kg 0.19 lbs	0.3 kgvol <sup>l</sup> 0.6 lbsvol <sup>l</sup>



## BIO THERM

The BioTherm range of High Performance 2-8°C solutions are fully compliant to IATA, ADR and DOT (49 CFR) transport regulations. These systems are qualified to ensure sample integrity for up to 96 hours against 2-8°C temperature exposures. 2-8°C temperature control is recommended when the samples being tested will be in transit for greater than 24 hours.

- 2-8°C systems qualified against the Intelsius High Performance Quad Profiles providing robust thermal protection all year round all required labelling for UN3373 Category B compliant shipping
- Sample integrity for up to 96 hours

### HIGH PERFORMANCE 2-8 SOLUTIONS

Product Code	Product Name	External Dimensions	Payload Dimensions	Mass	Volumetric Weight L/6	PathoShield 6/8 Capacity	PathoPak 1L Capacity
BT054P	BioTherm 15 HP 2-8	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	225 x 145 x 160 mm 8.8 x 5.7 x 6.3 in	6.0 kg 13.3 lbs	11.1 kgvol <sup>l</sup> 24.3 lbsvol <sup>l</sup>	2 / 7	0
BT061P	BioTherm 30 HP 2-8 Winter	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	275 x 190 x 220 mm 10.8 x 7.4 x 8.6 in	8.5 kg 18.6 lbs	17.0 kgvol <sup>l</sup> 37.3 lbsvol <sup>l</sup>	3 / 10	2
BT108	BioTherm 7 24 hour	308 x 308 x 310 mm 12.2 x 12.2 x 12.2 in	90 x 170 mm 3.5 x 6.7 in	3.5 kg 7.61 lbs	5.6 kgvol <sup>l</sup> 12.4 lbsvol <sup>l</sup>	N/A	1 - Supplied With System



## COMPLIANCE AND ORDERING

All DGP Intelsius UN3373 solutions are supplied ready to ship COVID-19 samples for testing. All systems are supplied with

- Pre-printed outers with full compliance labelling for UN3373 Category B infectious substances
- Absorbent material to minimise public risk in case of spillage or leaks
- Extra compliance labeling and security seals provided where required

Placing an order for packaging solutions can be completed through e-mail or by using the Intelsius e-commerce store. Please call if you would like to discuss your requirements.



<https://intelsius.com/product-category/coronavirus/>

T: +44 1904 607 390

E: [compliance@intelsius.com](mailto:compliance@intelsius.com)