

ORCA RESPONSE

2-8°C, 15-25°C, -20°C



TEMPERATURE
CONTROLLED
PACKAGING



ORCA Response is a multi-use high-performance temperature-controlled packaging solution designed for the transportation of packed red blood cells, platelets and pharmaceutical products. Payloads are maintained between +2°C and +8°C, +5°C and +15°C or below -20°C for a minimum of 40 hours utilising Phase Change Material (PCM) and Vacuum Insulation Panel (VIP) technologies.

A rectangular payload space reduces wasted space and minimises volumetric weight to allow solutions to be easily carried by hand for personal medication transport, direct to patient (DTP) shipments as well sent through courier networks. ORCA response is qualified against ISTA 7D Summer and Winter temperature profiles with a minimum performance of 40 hours.

Designed with the end user in mind ORCA Response combines high performance technology and ease of use where rapid response packaging is required. This optimum performance and small volumetric weight make ORCA Response a clear choice for temperature sensitive transportation and emergency applications.

Features

- 1L and 2L payload options
- Easy to assemble
- Simple coolant inventory
- Universal pack out
- Simple conditioning and assembly
- 3 temperature ranges: 2 to 8°C, 15 to 25°C and -20°C

For more information please contact Intelsius UK compliance@intelsius.com



ORCA RESPONSE



Code	Product Name	External Dimensions	Payload Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Minimum Duration
ORS111	ORCA Response 1L 2°C- 8°C	347 x 260 x 176mm 13.7 x 10.2 x 6.9in	190 x 130 x 44mm 7.5 x 5.1 x 1.7in	1.0L	4.4kg 9.7lbs	2.7 kgvol ^c 6.0 lbsvol ^c	40+
ORS221	ORCA Response 2L 2°C- 8°C	347 x 260 x 223mm 13.7 x 10.2 x 8.8in	190 x 130 x 91mm 7.5 x 5.1 x 3.6in	2.2L	5.4kg 11.9lbs	3.4 kgvol ^c 7.4 lbsvol ^c	60+
ORS112	ORCA Response 1L 15°C- 25°C	347 x 260 x 176mm 13.7 x 10.2 x 6.9in	190 x 130 x 44mm 7.5 x 5.1 x 1.7in	1.0L	4.7kg 10.4lbs	2.7 kgvol ^c 6.0 lbsvol ^c	40+
ORS222	ORCA Response 2L 15°C- 25°C	347 x 260 x 223mm 13.7 x 10.2 x 8.8in	190 x 130 x 91mm 7.5 x 5.1 x 3.6in	2.2L	5.8kg 12.8lbs	3.4 kgvol ^c 7.4 lbsvol ^c	60+
ORS113	ORCA Response 1L -20°C	347 x 260 x 176mm 13.7 x 10.2 x 6.9in	190 x 130 x 44mm 7.5 x 5.1 x 1.7in	1.0L	4.9kg 10.8lbs	2.7 kgvol ^c 6.0 lbsvol ^c	40+
ORS223	ORCA Response 2L -20°C	347 x 260 x 223mm 13.7 x 10.2 x 8.8in	190 x 130 x 91mm 7.5 x 5.1 x 3.6in	2.2L	6.0kg 13.2lbs	3.4 kgvol ^c 7.4 lbsvol ^c	40+

*Specifications subject to change

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8, Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0) 26 45699
salesireland@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai. Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com