

# BIOCARRIER 650

2-8°C, 15-25°C



**SAMPLE  
TRANSPORT  
PACKAGING**



**The Intelsius BioCarrier range is designed for local courier routes and multiple site sample collection or drug protection.**

When used with 95kPa primaries, the BioCarrier 650 acts as a leak proof secondary and outer packaging in order to meet regulatory requirements. Intelsius recommends the use of 95kPa compliant pouches (available separately) if your primaries are not 95kPa.

## THIS PRODUCT CAN BE USED FOR

- Multisite sample collections
- Diagnostic samples
- Screening programmes
- Vaccination delivery

## FEATURES

- Qualified for both 2-8°C and 15-25°C this product offers thermal protection for both temperature ranges with the same bill of materials
- Supplied pre labelled for Category B, Biological Substances
- Rigid reusable transit case with leak proof EPP liner

## BIOCARRIER 650 RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions	Payload Volume	Mass	Volumetric weight L/6	2-8°C Duration (hours)	15-25°C Duration (hours)
BC022	BioCarrier 4 650 8.0L	397 x 295 x 326 mm 15.6 x 11.6 x 12.8 in	200 x 200 x 200 mm 7.8 x 7.8 x 7.8 in	8.0L	4.9 kg 10.8 lbs	6.36 kgvol <sup>c</sup> 14.0 lbsvol <sup>c</sup>	15:30	24:00
BC026	BioCarrier 5 650 10.8L	397 x 295 x 425 mm 15.6 x 11.6 x 16.7	200 x 200 x 270 mm 7.8 x 7.8 x 10.6 in	10.8L	6.6 kg 14.6 lbs	8.3 kgvol <sup>c</sup> 18.3 lbsvol <sup>c</sup>	19:50	24:00

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



# BIOCARRIER 650

2-8°C & 15-25°C SOLUTIONS



**SAMPLE**  
TRANSPORT  
PACKAGING

**UNITED KINGDOM:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**GERMANY:** DGP Intelsius GmbH  
Eulerweg 11  
64291 Darmstadt  
Germany  
t: +49 6151 4936650  
f: +49 (0) 6151 4936651  
info@intelsius.de

**AMERICA:** DGP Intelsius LLC  
1399 Perry Rd, Suite 190,  
Plainfield, IN 46168,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**IRELAND:** DGP Intelsius Ireland Ltd.  
Unit 8. Coolcaslagh Industrial Estate,  
Coolcaslagh, Killarney  
Co. Kerry, V93 K006, Ireland  
t: +353 (0)64 668 61 71  
salesireland@intelsius.com

**MALAYSIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
Sun Srishti, C 10, Saki Vihar Road,  
Chandivali Powai. Mumbai - 400072  
Maharashtra - India  
+91 22 49746536  
compliance@intelsius.com