

BIOMAILER



**SAMPLE
TRANSPORT
PACKAGING**



The Intelsius BioMailer is a Category B compliant transport solution designed for use within national and global postal networks to ship small samples such as blood vials and slides.

THIS PRODUCT CAN BE USED FOR

- Category B, Biological Substances
- Screening programmes
- Clinical trials
- Swabs
- Vials

FEATURES

- Reusable rigid plastic case ensures samples are protected during transit
- Refurbishment kits available to replace single use components

BIOMAILER RANGE

Product Code	Product Name	External Dimensions	Internal Pouch Dimensions	Mass	Volumetric Weight L/6	Case Quantity
BM003	BioMailer with A5 Pouch	215 x 165 x 30 mm 8.5 x 6.5 x 1.2 in	155 x 180mm 6.1 x 7.0 in	0.2 kg 0.44 lbs	0.2 kgvol ^c 0.44 lbsvol ^c	50
BM004	BioMailer A5 Pouch Refurbishment Kit	N/A	155 x 180mm 6.1 x 7.0 in	0.04 kg 0.1 lbs	N/A	100
BM001	BioMailer with A5 Bag	215 x 165 x 30 mm 8.5 x 6.5 x 1.2 in	150 x 190 mm 5.9 x 7.4 in	0.2 kg 0.44 lbs	0.2 kgvol ^c 0.44 lbsvol ^c	50
BM002	BioMailer A5 Bag Refurbishment Kit	N/A	150 x 190 mm 5.9 x 7.4 in	0.04 kg 0.1 lbs	N/A	100

For more information please contact Intelsius UK cs@intelsius.com

BIOMAILER



SAMPLE
TRANSPORT
PACKAGING

UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8. Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0)64 668 61 71
salesireland@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai. Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com