

BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



The BioTherm range of Category A and Category B solutions come complete with the Intelsius PathoPak sample transport systems and are fully compliant to IATA, ADR and DOT (49 CFR) transport regulations. The systems are qualified to ensure sample integrity for up to 96 hours against 2-8°C, 15-25°C and dry ice. Each system is provided with handles and the correct compliance labelling to enable careful handling throughout the cold chain.

THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Complete temperature controlled Category A and Category B compliant overpack systems
- Range of systems to transport samples at 2-8°C, 15-25°C and deep frozen to maintain sample integrity
- 4GU compliant packaging systems for use with any primary receptacle
- Ergonomically engineered, reusable, sample bottle for safe manual handling and easy extraction

For more information please contact Intelsius UK cs@intelsius.com

BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE
TRANSPORT
PACKAGING**



DRY ICE SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Empty Mass	Dry Ice Mass	Volumetric Weight L/6	Duration
BT090	BioTherm 10A Dry Ice	357 x 309 x 315 mm 14.1 x 12.2 x 12.5 in	90 dia x 170 mm 3.5 dia x 6.7 in	1 L	1.4 kg 3.1 lbs	6.9 kg 15.3 lbs	5.8 kgvol ^c 12.8 lbsvol ^c	72
BT052	BioTherm 15A Dry Ice	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	2.4 kg 5.3 lbs	11.5 kg 25.4 lbs	11.1 kgvol ^c 24.4 lbsvol ^c	96
BT062	BioTherm 30A Dry Ice	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 255 mm 4.9 dia x 10.0 in	3 L	3.3 kg 7.1 lbs	25.0 kg 55.2 lbs	17.0 kgvol ^c 37.4 lbsvol ^c	96 +

HIGH PERFORMANCE 2-8°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT109	BioTherm 15 A 2-8 96 Hour	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	90 dia x 134 mm 3.5 dia x 5.3 in	800mL	6.19 kg 13.71 lbs	11.1 kgvol ^c 24.3 lbsvol ^c	96
BT056P	BioTherm 30A HP 2-8 96 Summer	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	10.2 kg 22.4 lbs	17 kgvol ^c 37.4 lbsvol ^c	96
BT057P	BioTherm 30A HP 2-8 96 Winter	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	9.0 kg 19.8 lbs	17 kgvol ^c 37.4 lbsvol ^c	96

CRT 15-25°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT110	BioTherm 15 A CRT 15-25	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	90 dia x 134 mm 3.5 dia x 5.3 in	800mL	7.79 kg 17.2 lbs	11.1 kgvol ^c 24.3 lbsvol ^c	96
BT105	BioTherm 30A CRT 15-25	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	11.9 kg 26.1 lbs	17 kgvol ^c 37.4 lbsvol ^c	96

UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8, Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0)64 668 61 71
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai, Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com