

ORCA SAMPLE TRANSPORT

2-8°C, 15-25°C, -20°C,

UN3373, UN2900, UN2814



TEMPERATURE
CONTROLLED
SAMPLE
PACKAGING



ORCA Sample Transport combines ORCA Single-Use, Multi-Use and Connect temperature-controlled systems with PathoPak sample transport packaging to give you a complete shipping solution for temperature sensitive samples.

Combining our PathoPak sample transport packaging and ORCA temperature-controlled systems into one, complete shipping solution for samples requiring temperature control provides a fully UN2900, UN2814 Category A Infectious Materials, and UN3373 Biological Substances compliant sample transport solution that ensures temperature stability throughout transit.

Features

- Fully compliant sample transport packaging solution
- High performance temperature-controlled overpack qualified against ISTA 7D temperature profiles
- Provided with all regulatory markings and labels to help ensure your dangerous goods compliance
- Available in a wide range of sizes
- Multi and single use solutions available
- Logging and tracking options available

For more information please contact Intelsius UK compliance@intelsius.com



CHOOSE YOUR ORCA

All ORCA overpack systems are pre-qualified and offer outstanding performance for a minimum 96 hours against ISTA 7D testing standards, meaning costly temperature excursions are avoided. When used as an overpack for infectious samples ORCA solutions should be labelled accordingly, a multi-label sheet (LA019) can be purchased separately.

ORCA Size Payload Volume	ORCA Single - Use		ORCA Multi - Use		ORCA Connect	
	Product Code	Temperature Range	Product Code	Temperature Range	Product Code	Temperature Range
ORCA 4L	OSU111	2-8°C	OMU111	2-8°C	OMXU111	2-8°C
	OSU112	15-25°C	OMU112	15-25°C	OMXU112	15-25°C
	OSU113	-15--25°C	OMU113	-15--25°C	OMXU113	-15--25°C
ORCA 11L	OSU221	2-8°C	OMU221	2-8°C	OMX221	2-8°C
	OSU222	15-25°C	OMU222	15-25°C	OMX222	15-25°C
	OSU223	-15--25°C	OMU223	-15--25°C	OMX223	-15--25°C
ORCA 23L	OSU331	2-8°C	OMU331	2-8°C	OMX331	2-8°C
	OSU332	15-25°C	OMU332	15-25°C	OMX332	15-25°C
	OSU333	-15--25°C	OMU333	-15--25°C	OMX333	-15--25°C
ORCA 50L	OSU441	2-8°C	OMU441	2-8°C	OMX441	2-8°C
	OSU442	15-25°C	OMU442	15-25°C	OMX442	15-25°C
	OSU443	-15--25°C	OMU443	-15--25°C	OMX443	-15--25°C

PATHOPAK FIT GUIDE

A strong, durable and highly protective secondary, when combined with ORCA overpack systems, PathoPak sample transport solutions provide the ideal secondary packaging for transporting temperature-sensitive samples such as blood, urine, tissue or swab samples, and come pre-printed with all the relevant UN markings. The table below shows how many of each size PathoPak fit inside your chosen ORCA solution.

PathoPak Size Payload Volume	Product Code	Case Quantity	ORCA Single-Use Multi-Use Connect 4L	ORCA Single-Use Multi-Use Connect 11L	ORCA Single-Use Multi-Use Connect 23L	ORCA Single-Use Multi-Use Connect 50L
PathoPak 800mL	PP034	20	1	4	4	18
PathoPak 1L	PP011	20	0	4	4	9
PathoPak 2L	PP012	10	0	1	1	4
PathoPak 3L	PP013	10	0	0	1	4

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8, Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0) 26 45699
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai, Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com