



**PathoPak sample transport packaging solutions are strong, durable, reusable, rigid secondary containers compliant with IATA, ADR, and CFR 49 (DOT) transport regulations for UN Class 6.2 Infectious Substances Category A and UN3373 Category B Infectious Substances.**

PathoPak packaging is suitable for shipping a wide range of primaries, including blood tubes, specimen containers, swabs, blood collection/transfusion bags and medical devices.

These reusable rigid systems are cost-effective and are available in four sizes: 800mL, 1L, 2L and 3L, allowing you to choose the right size for your payload, minimising wasted space.

## THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

## FEATURES:

- Complete Category A and Category B solutions include 95kPa tested secondary container, absorbent, bubble wrap, labels, and UN markings
- Reusable
- Ergonomically engineered bottle design for safe manual handling
- 4GU certified
- Suitable for cleaning by autoclave
- Full refurbishment services available

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)



**SAMPLE  
TRANSPORT  
PACKAGING**



## PATHOPAK RANGE

Product Code	Product Name	External Dimensions	Payload Dimensions (Bottle opening diameter x internal height)	Payload Volume	Mass	Volumetric Weight L/6
PP034	DGP PathoPak 800mL	112 x 110 x 150 mm 4.4 x 4.3 x 5.9 in	90 x 134 mm 3.5 x 5.3 in	800mL	0.19 kg 0.41 lbs	0.31 kgvol <sup>c</sup> 0.68 lbsvol <sup>c</sup>
PP011	DGP PathoPak 1L	112 x 110 x 190 mm 4.4 x 4.3 x 7.5 in	90 x 170 mm 3.5 x 6.7 in	1 L	0.24 kg 0.53 lbs	0.39 kgvol <sup>c</sup> 0.86 lbsvol <sup>c</sup>
PP012	DGP PathoPak 2L	153 x 156 x 193 mm 6.0 x 6.1 x 7.6 in	125 x 170 mm 4.9 x 6.7 in	2 L	0.45 kg 0.99 lbs	0.77 kgvol <sup>c</sup> 1.69 lbsvol <sup>c</sup>
PP013	DGP PathoPak 3L	160 x 160 x 280 mm 6.3 x 6.3 x 11.0 in	125 x 255 mm 4.9 x 10.0 in	3 L	0.6 kg 1.32 lbs	1.19 kgvol <sup>c</sup> 2.63 lbsvol <sup>c</sup>

**UNITED KINGDOM:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom

t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**AMERICA:** DGP Intelsius LLC  
1399 Perry Rd, Suite 190,  
Plainfield, IN 46168,  
United States of America

t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**MALAYSIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia

t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com

**GERMANY:** DGP Intelsius GmbH  
Eulerweg 11  
64291 Darmstadt  
Germany

t: +49 6151 4936650  
f: +49 (0) 6151 4936651  
info@intelsius.de

**IRELAND:** DGP Intelsius Ireland Ltd.  
Unit 8, Coolcaslagh Industrial Estate,  
Coolcaslagh, Killarney  
Co. Kerry, V93 K006, Ireland

t: +353 (0)64 668 61 71  
salesireland@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
Sun Srishti, C 10, Saki Vihar Road,  
Chandivali Powai, Mumbai - 400072  
Maharashtra - India

+91 22 49746536  
compliance@intelsius.com