

PATHOSEAL

SAMPLE TRANSPORT BAG



COMPONENT PACKAGING



Intelsius offer a range of flexible secondary packaging ideal for the transportation of a variety of primary containers such as blood vials, swabs, blood collection bags and medical devices.

The PathoSeal 95 range is designed and tested to withstand a pressure difference of 95 kPa to meet UN3373 Biological Substance, Category B regulations.

The PathoSeal Liquitite range is designed and tested to meet Category B regulations when used in conjunction with 95 kPa primaries.

All PathoSeal products feature a simple closure system enabling ease of packing, plus the unique tamper evident seal ensures the integrity of samples. In addition the bags are sequentially numbered, providing a chain of custody and additional security. An external document pouch is also present for all sizes with the exception of PathoSeal 95 A6.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Category B, Biological Substances
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Suitable for use within Category B solutions
- Tamper evident seal
- Sequentially numbered for sample traceability

For more information please contact Intelsius UK cs@intelsius.com



intelsius.com

PATHOSEAL

SAMPLE TRANSPORT BAG



COMPONENT
PACKAGING



PATHOSEAL 95 RANGE

Product Code	Product Name	Internal Dimensions (L x W)
BZ014	PathoSeal 95 A3	394 x 310 mm 15.5 x 12.2 in
BZ015	PathoSeal 95 A3 with 300ml Absorbent	394 x 310 mm 15.5 x 12.2 in
BZ006	PathoSeal 95 A4	284 x 200 mm 11.2 x 7.9 in
BZ009	PathoSeal 95 A4 with 200ml Absorbent	284 x 200 mm 11.2 x 7.9 in
BZ007	PathoSeal 95 A5	201 x 152 mm 7.9 x 6.0 in
BZ008	PathoSeal 95 A5 with 50ml Absorbent	201 x 152 mm 7.9 x 6.0 in
BZ020	PathoSeal 95 A6	139 x 99 mm 5.5 x 3.9 in
BZ013	PathoSeal 95 A6 with 50ml Absorbent	139 x 99 mm 5.5 x 3.9 in

PATHOSEAL LIQUITITE

Product Code	Product Name	Internal Dimensions (L x W)
BZ001	PathoSeal Liquitite 29/18	720 x 450 mm 28.3 x 17.7 in

UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8, Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0)64 668 61 71
salesireland@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai, Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com