

PATHOSHIELD



**SAMPLE
TRANSPORT
PACKAGING**



The PathoShield range are complete shipping solutions designed to carry a wide range of primaries, including blood vials, swabs and sample containers. Each system is pre-printed with the required markings and includes the relevant components to ship samples in compliance with transport regulations.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

FEATURES

- Complete Category A and Category B solutions include flexible 95kPa tested compliant secondary container, absorbent, labels, and UN markings
- Range of sizes suitable for shipping through postal and courier networks
- Carry a wide range of primaries including blood vials, swabs and sample containers
- Can be supplied as an insulated pack
- Supplied flat packed for efficient storage and shipping prior to use

For more information please contact Intelsius UK cs@intelsius.com



PATHOSHIELD RANGE

Product Code	Product Name	External Dimensions	Pouch Internal Dimensions (W x L)	Number of Pouches	Number of Tubes	Mass	Volumetric Weight L/6
Category A Solutions							
PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	308 x 230 mm 12.1 x 9.1 in	1	3	0.23 kg 0.5 lbs	0.9 kgvol ⁶ 2.1 lbsvol ⁶
Category B Solutions							
PA003	PathoShield 3 w A4 Pouch	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	308 x 230 mm 12.1 x 9.1 in	1	6	0.16 kg 0.35 lbs	0.5 kgvol ⁶ 1.0 lbsvol ⁶
PA002	PathoShield 4 w A4 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	308 x 230 mm 12.1 x 9.1 in	1	6	0.23 kg 0.5 lbs	0.9 kgvol ⁶ 2.1 lbsvol ⁶
PA001	PathoShield 5 w A5 Pouch	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.09 kg 0.20 lbs	0.4 kgvol ⁶ 0.8 lbsvol ⁶
PA006	PathoShield 6 w A5 Pouch	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	155 x 180mm 6.1 x 7.0 in	1	4	0.07 kg 0.15 lbs	0.2 kgvol ⁶ 0.4 lbsvol ⁶
PA007	PathoShield 7 w A5 Pouch	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	155 x 180mm 6.1 x 7.0 in	1	4	0.10 kg 0.21 lbs	0.3 kgvol ⁶ 0.6 lbsvol ⁶
PA008	PathoShield 8 w A5 Pouch	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	155 x 180mm 6.1 x 7.0 in	1	4	0.04 kg 0.09 lbs	0.1 kgvol ⁶ 0.2 lbsvol ⁶
PA023	PathoShield 3 w A4 Bag	191 x 145 x 102 mm 7.5 x 5.6 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.13 kg 0.28 lbs	0.5 kgvol ⁶ 1.0 lbsvol ⁶
PA024	PathoShield 4 w A4 Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	203 x 270 mm 7.9 x 10.6 in	1	-	0.19 kg 0.43 lbs	0.9 kgvol ⁶ 2.1 lbsvol ⁶
PA025	PathoShield 5 w A5 Bag	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.08 kg 0.18 lbs	0.4 kgvol ⁶ 0.8 lbsvol ⁶
PA026	PathoShield 6 w A5 Bag	191 x 145 x 43 mm 7.5 x 5.7 x 1.7 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.06 kg 0.13 lbs	0.2 kgvol ⁶ 0.4 lbsvol ⁶
PA027	PathoShield 7 w A5 Bag	191 x 145 x 57 mm 7.5 x 5.7 x 2.2 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.09 kg 0.19 lbs	0.3 kgvol ⁶ 0.6 lbsvol ⁶
PA028	PathoShield 8 w A5 Bag	191 x 145 x 22 mm 7.5 x 5.7 x 0.9 in	150 x 190 mm 5.9 x 7.4 in	1	-	0.04 kg 0.09 lbs	0.1 kgvol ⁶ 0.2 lbsvol ⁶
Outer							
PA014	PathoShield 4 Outer Box	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	-	-	-	0.16 kg 0.4 lbs	0.9 kgvol ⁶ 2.1 lbsvol ⁶
PA015	PathoShield 5 Outer Box	275 x 179 x 43 mm 10.8 x 7.0 x 1.7 in	-	-	-	0.07 kg 0.15 lbs	0.4 kgvol ⁶ 0.8 lbsvol ⁶
Mailing Bag							
MBUN3373C5*	Mailing bag UN3373 C5 1000 Units	240 x 170 mm 9.4 x 6.6 in	-	-	-	-	-

*UN3373 Mailing Bag is designed to be used in conjunction with a fully compliant UN3373 Category B packaging solution





PATHOSHIELD 4 INSULATED

Product Code	Product Name	External Dimensions	Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6
BT021	PathoShield 4 Insulated	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	251 x 151 x 62 mm 9.8 x 5.9 x 2.4 in	2.3 L	0.24 kg 0.53 lbs	0.9 kgvol ^c 2.1 lbsvol ^c

Product Code	Product Name	External Dimensions	Internal Dimensions	Number of Pouches	Number of Tubes	Mass	Volumetric Weight L/6
PA037	PathoShield 4 Insulated w A5 Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	150 x 190 mm 5.9 x 7.4 in	1	6	0.26 kg 0.57 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA038	PathoShield 4 Insulated w A5 Bag & Gel Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	150 x 190 mm 5.9 x 7.4 in	1	6	0.66 kg 1.45 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA039	PathoShield 4 Insulated w A5 Pouch	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	155 x 180 mm 6.1 x 7.0 in	1	6	0.26 kg 0.57 lbs	0.9 kgvol ^c 2.1 lbsvol ^c
PA040	PathoShield 4 Insulated w A5 Pouch & Gel Bag	291 x 191 x 102 mm 11.4 x 7.5 x 4.0 in	155 x 180 mm 6.1 x 7.0 in	1	6	0.66 kg 1.45 lbs	0.9 kgvol ^c 2.1 lbsvol ^c

UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom

t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America

t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia

t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany

t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8. Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland

t: +353 (0)64 668 61 71
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai. Mumbai - 400072
Maharashtra - India

+91 22 49746536
compliance@intelsius.com