

PATHOSORB ABSORBENT



COMPONENT
PACKAGING



The Intelsius PathoSorb is a Super Absorbent Polymer (SAP) that ensures excellent absorption capacity to size ratio. Upon absorption during accidental spillage, PathoSorb acquires a 'gel like' consistency, locking in the spilled material for easy disposal. PathoSorb is lint free, reducing sterile contamination risk.

Category A and Category B transport regulations specify that absorbent material must be placed between primary and secondary containers in sufficient quantity to absorb the entire contents of the primary containers.

FEATURES

- ▶ PathoSorb utilises SAP technology to provide effective absorbency whilst maximising available payload space.
- ▶ For multiple sample transportation, PathoSorb 4 & 6 Bay pouches provide a dual purpose; simultaneously protecting and separating samples.

For more information please contact Intelsius UK cs@intelsius.com



PATHOSORB ABSORBENT



COMPONENT
PACKAGING



PATHOSORB ABSORBENT RANGE

Product Code	Product Name	Case Quantity	Dimension	Water Absorbent
AB001	PathoSorb 50ml	1000	150 x 75mm	56
AB010	PathoSorb 100ml	1000	160 x 140mm	112
AB009	PathoSorb 200ml	100	225 x 175mm	223
AB002	PathoSorb 300ml	1000	305 x 220mm	335
AB007	PathoSorb 600ml	100	470 x 295mm	610
AB005	PathoSorb 4 Bay Pouch	1	225 x 125mm	155
AB008	PathoSorb 6 Bay Pouch	100	325 x 125mm	235

EUROPE: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICAS: DGP Intelsius LLC
7696 Zionsville Road,
Indianapolis, IN 46268,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

ASIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com



IRELAND: DGP Intelsius Teo
Unit 8, Ballyvourney Industrial Estate,
Ballyvourney, Co. Cork,
P12 A527, Ireland
t: +353 (0) 26 456 99
f: +353 (0) 26 456 98
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
2, Palak Pluto, Off Link Road,
Kanderpada – Dahisar (West),
Mumbai 400 068, Maharashtra, India
t: +91 22 2892 3675
f: +91 22 2892 3675
compliance@intelsius.com