



The Intelsius BioShield Range of unique patented systems are designed to meet all handling processes from sample extraction to laboratory presentation. BioShields are highly protective single-use shipping solutions. They are available for a variety of sample capacities, helping to reduce waste and cut costs associated with sample transportation.

The systems are well suited to a wide variety of Category B samples, including solids, liquids, and samples held in fixative. The BioShield helps reduce total costs by ensuring that all samples will arrive at the same time for a simple and efficient sample reception.

THIS PRODUCT CAN BE USED FOR

- Diagnostic samples
- Screening programmes
- Eradication schemes

FEATURES

- Complete Category B compliant solutions including 95kPa 90ml sample pots, tray, extraction utensils and leak proof bag
- Techni-Ice gel sheet available with next day systems for temperature stability during transport

For more information please contact Intelsius UK cs@intelsius.com



BIOSHIELD 'SAME DAY' RANGE

Product	Product Code (Inc. Pots)	Product Code (Exc. Pots)	Product Code (Inc. Pots, Spoon & Forceps)	Sample Capacity	External Dimensions	Volumetric Weight
BioShield 6	BSE061	BSE072	BSE055	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield 10	BSE075	BSE077	BSE038	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield SR 20	BSE056	BSE065	BSE057	20	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield SR 40	BSE058	BSE066	BSE059	40	519 x 340 x 125 mm 20.5 x 13.4 x 5 in	3.7 kgvol ^c 8.1 lbsvol ^c

BIOSHIELD 'NEXT DAY' RANGE

Product	Product Code (Inc. Pots)	Product Code (Exc. Pots)	Product Code (Inc. Pots, Spoon & Forceps)	Sample Capacity	External Dimensions	Volumetric Weight
BioShield 6	BSE062	BSE071	BSE051	6	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield 10	BSE074	BSE076	BSE073	10	291 x 191 x 102 mm 11.5 x 7.6 x 4.1 in	0.94 kgvol ^c 2.1 lbsvol ^c
BioShield SR 20	BSE050	BSE063	BSE064	20	524 x 340 x 107 mm 20.6 x 13.4 x 4.2 in	3.7 kgvol ^c 8.1 lbsvol ^c
BioShield SR 40	BSE041	BSE042	BSE067	40	524 x 340 x 107 mm 20.6 x 13.4 x 4.2 in	3.7 kgvol ^c 8.1 lbsvol ^c

SAMPLE EXTRACTION TOOLS

Product	Product Code	Dimensions	Quantity
Curved Bovine 3 Bar Spoon	SPO01	191mm x 29mm x 25mm	1000
Red Ovine Spoon	SPO03	224.9mm x 16.7mm x 8.2mm	1000
Cerebellum Ovine Atypical 2 Bar Spoon	SPO10	274.7mm x 16.1mm x 28.6mm	1000
Forceps Wide Tip Green	SPO12	130mm x 30mm x 13mm	20



UNITED KINGDOM: DGP Intelsius Ltd
1 Harrier Court, Airfield Business Park,
Elvington, York, YO41 4EA,
United Kingdom
t: +44 (0) 1904 607 390
f: +44 (0) 1904 607 391
cs@intelsius.com

AMERICA: DGP Intelsius LLC
1399 Perry Rd, Suite 190,
Plainfield, IN 46168,
United States of America
t: +1 317 452 4006
f: +1 317 290 4621
americasorders@intelsius.com

MALAYSIA: DGP Intelsius Sdn Bhd
Plot 11, Lorong Bemban 2, Estet
Perindustrian Bemban, 31000
Batu Gajah, Malaysia
t: +605 3666 399
f: +605 3669 399
compliance@intelsius.com

GERMANY: DGP Intelsius GmbH
Eulerweg 11
64291 Darmstadt
Germany
t: +49 6151 4936650
f: +49 (0) 6151 4936651
info@intelsius.de

IRELAND: DGP Intelsius Ireland Ltd.
Unit 8. Coolcaslagh Industrial Estate,
Coolcaslagh, Killarney
Co. Kerry, V93 K006, Ireland
t: +353 (0)64 668 61 71
salesireland@intelsius.com

INDIA: DGP Intelsius Cold Chain Pvt. Ltd.
Sun Srishti, C 10, Saki Vihar Road,
Chandivali Powai. Mumbai - 400072
Maharashtra - India
+91 22 49746536
compliance@intelsius.com